LifePhotonic

Instructions for Use (IFU)

Cylindrical Diffusor CD-2S3 / CD-403 / CD-4L3

Front Diffusor FD-4

Intended Purpose and Indications

The LifePhotonic Cylindrical and Front Diffusers are sterile, single use optical fibers intended to deliver spatially homogeneous laser light to biological tissue during minimally invasive endoscopic or intra operative procedures. By enabling predictable dosimetry, the Diffusers support photodynamic, photothermal or fluorescence guided interventions in

- Lung
- Oesophagus
- Head and Neck
- Bile Duct
- Prostate
- Urothelial tract
- Pancreas

with the Front Diffuser being limited to superficial, non-contact illumination in Head and Neck.

The devices themselves provide **no pharmacological, immunological or metabolic effect**, and overall treatment outcome depends among others on photosensitizer uptake, wavelength selection, total fluence and tissue oxygenation, which lie outside the scope of the Diffuser product. Before use, users must read and follow the SmPC and the treatment/study protocol of all photosensitive medicines involved in the treatment as well as the IFU of the laser and of the catheter or endoscope.

The intended patient population is limited to adults.

Contraindications

The risk class of the devices does not allow an application in contact with central circulatory system and with central nervous system.

Safety Instructions and warnings

LP_IFU-Diff-EDSFI_R; Rev. R, 2026-02-20

The diffusor may only be used by physicians trained in phototherapies.

The diffusor may only be used with a photosensitive medicine, with a laser, and a catheter or endoscopic insertion device, each either approved or used under study conditions in the targeted indication. The user is responsible for verifying whether the laser is suitable to work with the diffusor and whether the diffusor can be inserted into the catheter/endoscope. Read and follow the IFU of the laser system, and the catheter / endoscope. Read and follow the SmPC and the treatment/study protocol of the photosensitive medicine before use.

 Sterilized with EO. The EO indicator must be yellow. Adhere to regulations for handling sterile products.

 Single sterile barrier in outer pouch

 Do not use when packaging is opened or damaged. Do not use this product if you have doubts about its sterility.

 Single-Use  Medical Device

 Adhere to regulations for handling laser radiation. Ensure that all users are aware of handling fiber optics. Always wear protective eyewear during laser light delivery. Avoid eye or skin exposure to direct or scattered laser radiation. Provide eye protection to patient and personnel in the laser area.

 Users must be trained in handling diffusors. Contact LifePhotonic in case of training needs.

 Fragile device. The optical fiber may break if handled roughly or bent sharply! Pay attention to the smallest permissible bending diameter. Do not advance the diffusor against resistance.

 The diffusor must not contact metallic surfaces when emitting light to avoid potentially dangerous scattered laser radiation or a damage of the diffusor.

 Keep the optical connector clean and dry. Avoid touching the connector surface as it may stain or scratch the optical surface of the connector.

 Do not exceed maximum output power specifications.

Users should report any serious incident that has occurred in relation to the device to the manufacturer and to the competent authority of the EU member state in which the user and/or patient is established. In case the device is used in clinical studies or clinical investigations,

investigators shall report any serious adverse events and potential device deficiencies to the sponsor and to the authorities as applicable following study protocol instructions and local regulations. The sponsor is advised to report serious adverse events with causal or potential causal relationship with the device and potential device deficiencies to the manufacturer in a timely manner, as applicable by local regulations.

Device Characteristics Cylindrical Diffusor

| CD- | 2S3-xx-z | 403-xx-z | 4L3-xx-z |
|-----------------------------|------------------|------------------|----------|
| MECHANICAL | | | |
| Tip diameter | 0.75 mm | 0.98 mm | < 1.2 mm |
| Overall length of fiber | 3.0 m | | |
| Short-term bending diameter | 25 mm | 40 mm | |
| OPTICAL | | | |
| Active diffusor length xx | xx is 05 – 40 mm | xx is 10 - 90 mm | |
| Uniformity of diffusor | ± 20% | | |
| Max. power density | 600 mW/cm | | |
| Max. power | 2.5 W | | |
| Wavelength | 600 nm – 800 nm | | |
| Numerical Aperture | 0.37 | | |
| Core diameter | 200 µm | 400 µm | |
| Fiber connector | SMA 905 | | |

Options for X-ray markers are: with 2 X-ray markers (z = -MX), with 1 X-ray-marker (z = -M), or without X-ray-marker (z = -0).

Device Characteristics Front Diffusor

| | FD-4 |
|-----------------------------|----------|
| MECHANICAL | |
| Tip diameter | < 750 µm |
| Overall length of fiber | 3.0 m |
| Short-term bending diameter | 40 mm |
| OPTICAL | |
| Typical transmission | > 85% |

| | |
|--------------------|-----------------|
| Max. power | 3.0 W |
| Wavelength | 600 nm – 800 nm |
| Numerical Aperture | 0.37 |
| Core diameter | 400 µm |
| Fiber connector | SMA 905 |

Instructions for proper use

- Use the diffusor only with lasers with a SMA-compatible fiber port, accepting fibers with the respective numerical aperture and core diameter or smaller.
- Remove the device from its packaging and inspect the fiber for visible signs of damage. Securely attach the optical connector to the laser fiber port. Inspect the uniformity of the diffusor using the aiming beam of the laser system. Do not use the diffusor in case of obvious damage or if the light output does not appear uniform.
- Test the compatibility of the laser device by connecting the diffusor to the laser and successfully completing the automatic calibration routine of the laser, following the procedure described in the laser system's instructions for use. If the laser does not provide an automatic calibration routine or the automatic calibration is aborted, the laser is not suitable.
- Treatments with the device CD-2S3-xx have to be performed only in liquid environment like water or cooling fluid in order to provide sufficient cooling.

In case of unexpected calibration results, e.g. lower transmission than expected, do not use the diffusor.

If required by the treatment or study protocol, perform the calibration procedure under sterile conditions, e.g. by using a sterile, transparent cuvette or catheter in the calibration unit.

- Follow the SmPC and the treatment or study protocol of the photosensitive medicine to determine the required light dose. Set and if required calibrate the laser power accordingly following the IFU of the laser system.

If the diffusor shall be inserted via a catheter or endoscope, it must be used and positioned in the organ to be treated with the catheter or endoscope described in the instructions of the treatment or study protocol.

Avoid tissue residue, blood and/or expectoration on the diffusor. This may cause excessive light absorption, and thus may alter the treatment outcome and may damage the diffusor

Avoid inadvertent photo activation in non-target tissue. If necessary shield surrounding tissue from laser-light.

Storage

Protect from sunlight.

Keep dry.

Disposal

The product must be disposed of in accordance with local regulations for contaminated materials. Defective contaminated material may only be returned after disinfection.

The following risks may arise in case of an impermissible reprocessing or re-sterilization:

Infection due to insufficient sterility and/or biological contamination. Reprocessing may lead to invisible damages of the device and thus insufficient functionality and unintended treatment outcome.

Disclaimer of Warranties

LifePhotonic disclaims any liability in the case of use of the device without compliance to the instructions for use or for a purpose other than the intended use indicated in these instructions.

LifePhotonic liability is limited to the replacement of defective pieces, and to the damage, that is a direct consequence of the use of defective devices.

In no case *LifePhotonic* can be held responsible for an incorrect light dose applied during a phototherapy treatment.

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LifePhotonic GmbH